禁忌症及注意事项详见说明书 苏械广审(文)第281227-04003号

技术参数

		话田空刺哭		
吻合器型号	臂杆长度	适用穿刺器		
BQS45/60/45A/60A	190mm	12		
BQM45/60/45A/60A	250mm	12		
BQL45/60/45A/60A	350mm	12		

		計杆长度 适用穿刺器		
吻合器型号	臂杆长度	适用穿刺器		
BCS45/60/45A/60A	190mm	12		
BCM45/60/45A/60A	250mm	12		
BCL45/60/45A/60A	350mm	12		

钉仓型号	钉仓颜色	钛钉高度	闭合钉高度	备注
BQJN45、BRJN45、BCJN45	灰	2.0mm	0.75mm	
BQJW45、BRJW45、BCJW45	白〇	2.5mm	1.0mm	
BQJB45、BRJB45、BCJB45	蓝	3.5mm	1.5mm	
BQJY45、BRJY45、BCJY45	金	3.8mm	1.8mm	
BQJG45、BRJG45、BCJG45	绿 ●	4.1mm	2.0mm	
BQJH45、BRJH45、BCJH45	黑 ●	4.4mm	2.3mm	
BQJN45K、BRJN45K、BCJN45K	灰	2.0mm	0.75mm	2/3
BQJW45K、BRJW45K、BCJW45K	白 O	2.5mm	1.0mm	Sils
BQJB45K、BRJB45K、BCJB45K	蓝	3.5mm	1.5mm	Selfs.
BQJY45K、BRJY45K、BCJY45K	金	3.8mm	1.8mm	2 ^N Z
BQJG45K、BRJG45K、BCJG45K	绿 •	4.1mm	2.0mm	292
BQJH45K、BRJH45K、BCJH45K	黑 ●	4.4mm	2.3mm	Syl
BQJN60、BRJN60、BCJN60	灰	2.0mm	0.75mm	
BQJW60、BRJW60、BCJW60	白 O	2.5mm	1.0mm	
BQJB60、BRJB60、BCJB60	蓝	3.5mm	1.5mm	
BQJY60、BRJY60、BCJY60	金	3.8mm	1.8mm	
BQJG60、BRJG60、BCJG60	绿 •	4.1mm	2.0mm	
BQJH60、BRJH60、BCJH60	黑	4.4mm	2.3mm	
BQJN60K、BRJN60K、BCJN60K	灰	2.0mm	0.75mm	SIB
BQJW60K、BRJW60K、BCJW60K	白〇	2.5mm	1.0mm	25
BQJB60K、BRJB60K、BCJB60K	蓝	3.5mm	1.5mm	SIB
BQJY60K、BRJY60K、BCJY60K	金	3.8mm	1.8mm	SIB
BQJG60K、BRJG60K、BCJG60K	绿 •	4.1mm	2.0mm	SIB
BQJH60K、BRJH60K、BCJH60K	黑 ●	4.4mm	2.3mm	292

邦尼医疗科技(常州)有限公司

Bonnie Medical Device Technology (Changzhou)Co.,Ltd.

网址: www.bnylkj.com 即编: bnylkj2018@163.com

■ 地址: 江苏省常州市武进区湖塘镇湖塘科技产业园A1号 电话: 0519−8818 8331





一次性电动腔镜用 Disposable Powered Endoscopic

直线型切割吻合器及钉仓组件

Linear Cutter Stapler and Reload

